


FESCA-200 Femtosecond streak camera C6138 series




The FESCA-200 is an ultra fast streak camera with a temporal resolution of 200 femtoseconds (typ.). It is designed for use with single-shot or slow-repetitive phenomena. It can analyze the process of energy relaxation and the dynamics of chemical reaction in the femtosecond region in combination with femtosecond pulse laser.

Features

- Single-shot 200 fs (typ.) temporal resolution
- Simultaneous measurement of light intensity on both the temporal and spatial axis
- Dedicated readout system



▲ Streak image and intensity profile of light pulses from Ti:Sapphire laser measured with the FESCA-200




▲ This streak image, taken with the FESCA-200, was used to measure Cherenkov radiation. LINAC was used to generate an electronic pulse, which, in contact with the air, produced the Cherenkov radiation. Photo courtesy of Dr. Mitsuru Uesaka, Associate Professor, Nuclear Engineering Research Laboratory, Faculty of Engineering, University of Tokyo.

Applications

- Research of the process of energy relaxation of quantum well semiconductors
- Research of the dynamics of chemical reactions in the femtosecond region
- Research of the dynamics of ultra fast laser diodes, ultra fast optical logic devices, etc.
- Diagnosis of femtosecond lasers

System configuration



Operating principle


The light pulse to be measured is projected onto the slit and is focused by the lens into an optical image on the photocathode of the streak tube. Changing the temporal and spatial offset slightly each time, four light pulses, each with a different light intensity, are introduced through the slit and conducted to the photocathode.

Here, the photons are converted into a number of electrons proportional to the intensity of the incident light. The four light pulses are converted sequentially to electrons which are then accelerated and conducted towards the photocathode.

As the group of electrons created from the four light pulses passes between a pair of sweep electrodes, a high-voltage is applied (see above), resulting in a high-speed sweep (the electrons are swept in the direction from top to bottom). The electrons are deflected at different times, and at slightly different angles in the perpendicular

direction, and are then conducted to the MCP (micro-channel plate). As the electrons pass the MCP, they are multiplied several thousands of times and are then bombarded against the phosphor screen, where they are converted back into light.

The fluorescence image corresponding to the first incident light pulse is positioned at the top of the phosphor screen, followed by the others, with images proceeding in descending order; in other words, the axis in the perpendicular direction on the phosphor screen serves as the temporal axis. The brightness of the various fluorescence images are proportional to the intensities of the corresponding incident light pulses. The positions in the horizontal direction on the phosphor screen correspond to the positions of the incident light in the horizontal direction.



Operating Principle of the Streak Tube

Specifications

① FESCA-200 Femtosecond streak camera C6138-11, -12

• Main unit

| | | |
|------------------------------------|---|-------------------------|
| Temporal resolution | 200 fs (better than 300 fs) typ. (at the fastest sweep range) | |
| Effective phosphor screen size | Time axis | Approx. 10 mm |
| | Spatial axis | Approx. 9 mm |
| Sweep time/full screen (10 mm) | 20 ps, 50 ps, 100 ps, 200 ps, 500 ps, 1 ns | |
| Trigger jitter | Less than ± 20 ps | |
| Trigger delay | Approx. 30 ns (at the fastest sweep range) | |
| Maximum sweep repetition frequency | 100 Hz | |
| Trigger input | Maximum input voltage | ± 5 V / 50 Ω |
| | Maximum input repetition frequency | 150 MHz |
| | Trigger level | -4 V to +4 V |
| Power supply | AC 100 V to AC 240 V, 50 Hz/60 Hz | |
| Operating temperature | 0 °C to +40 °C | |
| Operating humidity | Less than 70 % (with no condensation) | |

• Streak tube

| | | |
|----------------------------------|----------------------------|-------------------|
| Photocathode/Material of Windows | Multi-alkali / Kovar glass | |
| Effective photocathode size | Max. 3 mm (H) | |
| Spectral response | C6138-11 | 280 nm to 850 nm |
| | C6138-12 | 300 nm to 1050 nm |
| Phosphor screen size | ϕ 18 mm | |
| Phosphor screen | P-43 | |
| Image magnification | 1 : 3 | |

• Image intensifier

| | |
|--------------------------------|----------------|
| Effective photocathode size | ϕ 18 mm |
| Photocathode | GaAsP |
| Light gain | More than 2500 |
| Phosphor screen | P-43 |
| Effective phosphor screen size | ϕ 18 mm |

② Data analyzer C6760-36

• Data analysis module M5755-36

| | |
|-------------------------|---|
| Image Input board | IEEE1394b |
| I/F | PCI |
| Data acquisition | Monitoring, Analog integration, Photon counting |
| External device control | Streak camera, Spectrometer etc. |
| Profile analysis | Realtime monitoring, FWHM, Peak count etc. |
| Calibration | Dark current, Shading, Time axis, Wavelength axis, Jitter |
| File format (image) | Binary (16 bit), TIFF, ASCII |
| File format (profile) | TEXT |
| Data processing | Windows PC |

③ ORCA[®]-R2 Digital CCD camera C10600-10B

| | | |
|------------------------------|---|--------------------------------------|
| Effective number of pixels | 1344(H) \times 1024(V) | |
| Cell size | 6.45 μ m (H) \times 6.45 μ m (V) | |
| Frame rate | Normal scan | 8.5 Hz (binning : 1 \times 1) |
| | | 15.6 Hz (binning : 2 \times 2) |
| | Fast scan | 16.2 Hz (binning : 1 \times 1) |
| | | 28.4 Hz (binning : 2 \times 2) |
| Readout noise (rms) typ. | Normal scan | 6 electrons |
| | Fast scan | 10 electrons |
| Dynamic range typ. | 3000 : 1 (Normal scan / binning : 1 \times 1) | |
| Cooling method / temperature | Forced-air cooled | - 35 °C |
| | Water cooled | - 40 °C (Water temperature : +20 °C) |
| Dark current | 0.0005 electrons/pixel/s (- 40 °C) | |
| Dual A/D converter | 12 bit or 16 bit | |
| Exposure time | 10 μ s to 4200 s | |


④ Input optics

| | | |
|------------------------|-----------------------|---------------------|
| Type number | Input optics A1976-01 | Mirror optics A6856 |
| Spectral transmittance | 200 nm to 1600 nm | 200 nm to 1600 nm |
| Image magnification | 1 : 1 | 1 : 1 |
| Effective F number | 5.0 | 4.0 |


< Others >

- Pin diode head C1083-01
- Delay unit C1097-05

● C6138 main unit (Approx. 26 kg)



● C6138 power supply unit (Approx. 2.6 kg)



● The temporal resolution of FESCA-200

The temporal resolution of FESCA-200 can be limited depending on wavelength or input light intensity. (See Figure 1, 2)

The photoelectron from photocathode spread spatially along as the wavelength become shorter. Therefore, the temporal resolution is limited by wavelength.

The photoelectron can be spread spatially, each photoelectron bound back while traveling inside streak tube if a number of photoelectron is increased. Thus, the temporal resolution can be limited.




Figure 1




Figure 2

★ ORCA is registered a trademark of Hamamatsu Photonics K.K.

★ Product and software package names noted in this documentation are trademarks or registered trademarks of their respective manufacturers.

● Subject to local technical requirements and regulations, availability of products included in this promotional material may vary. Please consult your local sales representative.

● Information furnished by HAMAMATSU is believed to be reliable. However, no responsibility is assumed for possible inaccuracies or omissions.

Specifications and external appearance are subject to change without notice.

© 2015 Hamamatsu Photonics K.K.

HAMAMATSU PHOTONICS K.K. www.hamamatsu.com

HAMAMATSU PHOTONICS K.K., Systems Division

812 Joko-cho, Higashi-ku, Hamamatsu City, 431-3196, Japan, Telephone: (81)53-431-0124, Fax: (81)53-435-1574, E-mail: export@sys.hpk.co.jp

U.S.A.: Hamamatsu Corporation: 360 Foothill Road, Bridgewater, N.J 08807, U.S.A., Telephone: (1)908-231-0960, Fax: (1)908-231-1218 E-mail: usa@hamamatsu.com

Germany: Hamamatsu Photonics Deutschland GmbH.: Arzbergerstr. 10, D-82211 Herrsching am Ammersee, Germany, Telephone: (49)8152-375-0, Fax: (49)8152-265-8 E-mail: info@hamamatsu.de

France: Hamamatsu Photonics France S.A.R.L.: 19, Rue du Saule Trapu, Parc du Moulin de Massy, 91882 Massy Cedex, France, Telephone: (33)1 69 53 71 00, Fax: (33)1 69 53 71 10 E-mail: infos@hamamatsu.fr

United Kingdom: Hamamatsu Photonics UK Limited: 2 Howard Court, 10 Tewin Road, Welwyn Garden City, Hertfordshire AL7 1BW, UK, Telephone: (44)1707-294888, Fax: (44)1707-325777 E-mail: info@hamamatsu.co.uk

North Europe: Hamamatsu Photonics Norden AB: Torshamnsgatan 35 16440 Kista, Sweden, Telephone: (46)8-509-031-00, Fax: (46)8-509-031-01 E-mail: info@hamamatsu.se

Italy: Hamamatsu Photonics Italia S.r.l.: Strada della Moia, 1 int. 6 20020 Arese (Milano), Italy, Telephone: (39)02-93581733, Fax: (39)02-93581741 E-mail: info@hamamatsu.it

China: Hamamatsu Photonics (China) Co., Ltd.: B1201 Jiaming Center, No.27 Dongsanhuan Beilu, Chaoyang District, Beijing 100020, China, Telephone: (86)10-6586-6006, Fax: (86)10-6586-2866 E-mail: hpc@hamamatsu.com.cn