



the filler for every solution

Continual **physico-chemical studies** on cross-linking have resulted in the synthesis of all the experience gained to date, resulting in the creation and development of an **ingenious** product that incorporates all the benefits desired by medical staff and patients, **i.e. safe, natural, long-lasting and comfortable implants**, in a single solution: for this reason Genyal Polyvalent, the polyvalent filler for all types of cosmetic correction, was formulated.

Genyal® Polivalent:

- treats all types of superficial and medium-depth wrinkles
- remodels the lip contours and increases lip volume

SPECIAL FEATURES

The latest advances in the field have shown that polyvalent fillers are, currently, the best possible upgrade but Genyal's special features make it a unique product:

- **non-animal origin:** sodium hyaluronate synthesised by bio-fermentation dual-action, high-concentration filler: 23mg/1ml of cross-linked hyaluronic acid with bio-stimulating free hyaluronic acid
- **ultra high molecular weight product (3 million Daltons):** this characteristic, together with cross-linkage, allows us to increase the resistance to degradation by hyaluronidases to a maximum and thus maximise the duration of the product in the tissue.
- **cross-linking with butanediol diglycidyl ether (BDDE),** the cross-linking agent universally recognised as the least toxic of all, which ensures that the risk of allergies is minimal.
- **single-phase formula:** all the products in the **Genyal®** range are single-phase and **therefore contain no particles, which ensures the homogeneity of the gel.**

HYALURONIC DERMAL FILLER

GENYAL®

Polyvalent



CONTENT	1ML
CONCENTRATION HA	23MG/ML
RETICULATING AGENT	BDDE
DEGREE OF CROSSLINKING	★★★★
DINAMIC VISCOSITY	500.000 mPas
NEEDLES	2 X 27G

TREATMENTS: lip contour, barcode, glabella, facial wrinkles, forehead wrinkles, periorcular wrinkles, naso genian folds, perilabial wrinkles, labial commissures.





HA-GLYCEROL BIO-RESTRUCTURING GEL

GENYAL[®]

Genyalift



CONTENT	1ML
CONCENTRATION HA	18MG/ML
STABILISING AGENT	GLYCEROL 20MG/ML
MOLECULAR WEIGHT	3 MLN DALTON
NEEDLES	2 X 30G

TREATMENTS: glabella, forehead wrinkles, periocular wrinkles, cheeks, naso genian folds, perilabial wrinkles, labial commissures, neck, facial wrinkles, hands, stretch marks, décolleté, face.

bio-stimulation at the highest level

It is now indisputable that filling wrinkles and increasing the volume in the face is not the only way to obtain total rejuvenation of the face; this is why natural, non-cross-linked hyaluronic acid has an important role in combating and preventing skin aging, allowing skin elasticity and firmness to be regained.

The particular formulation of **Genyal® Genyalift**, with its **extremely high concentration of hyaluronic acid (18 mg/1 ml) together with a high concentration of glycerol (20 mg/ml) and ultra high molecular weight (3 million Daltons: double or even three times that of many rival formulas!)**, allows the recreation of the optimum physiological environment for fibroblast proliferation and migration and for the neo-synthesis of collagen, elastin and other constituents of the extracellular matrix.

Genyal® Genyalift can be used not only on all areas of the face, but also in other regions such as the neckline and the backs of the hands, with different programmes depending on the degree of skin aging; generally, 2-3 initial treatments are planned at short time intervals (15-30 days) and then other maintenance sessions are continued about 2-3 months apart.



for all types of volume



Genyal® Volumae is the filler with the best balance between high viscosity, elasticity and long lasting. Combined to total confort, satisfaction and reliability during filler implantation, for an excellent, natural, finished filling effect.

The excellent results achieved with **Genyal® Polyvalent** led to the creation of a specific product that gives volume to the face while maintaining the characteristics of the **Genyal®** range: the result is a highly cross-linked gel (also containing 23 mg/1 ml of hyaluronic acid) which is injected with a 27G needle to reach the desired level of filling without causing trauma to the tissues.

Genyal® Volumae:

- fills lines and treats deep wrinkles in thick skin
- restores or creates facial volume
- increases the volume of the lips

The duration is exceptional for a reabsorbable product: **it varies between 8 and 18 months following implantation.**

SPECTACULAR EFFECTS AND EXCEPTIONAL DURATION

HYALURONIC DERMAL FILLER

GENYAL®

Volumae



CONTENT	1ML
CONCENTRATION HA	23MG/ML
RETICULATING AGENT	BDDE
DEGREE OF CROSSLINKING	★★★★
DINAMIC VISCOSITY	1MLN mPas
NEEDLES	2 X 27G

TREATMENTS: cheeks, naso genian folds, labial commissures, increase in lip volume, chin, facial wrinkles, nasal correction.





Xcelens Exclusive Formula

Xcelens International is pleased to present an innovative and exclusive range of hyaluronic acid products, uniquely developed to guarantee medical staff, and therefore patients, total comfort, satisfaction and reliability during filler implantation and bio-stimulation treatments: **GENYAL®!**

The **GENYAL®** range is developed by an expert **Swiss-Italian team** of chemists, technicians, engineers and professionals operating in the field of aesthetic medicine with a view to creating more advanced and technologically avant-garde solutions year after year.

- **Monophasic**
- **Easy to handle**
- **Immediate results**
- **Natural effects**
- **Highly purified HA**
(Lowest level of endotoxin < 0.05 EU/g)
- **Great visco-elasticity balance**
- **Natural reabsorption**
- **Long-lasting**
(8 to 18 months)

Innovative technology



go to www.genyal.ch
to see our technology